



Be Right™



Bottle, m-ColiBlue24 broth, 100 mL (50 test)

Artikel-Nr.: 2608442

EUR Preis: Kontakt

Lieferung innerhalb von 2 Wochen

m-ColiBlue24 Broth in 100 mL Bottle

In 24 hours, Hach's m-ColiBlue24 Broth gives you simultaneous results for total coliforms and E. coli. Hach's ready-to-use media eliminates measuring, mixing, and autoclaving steps necessary to prepare dehydrated media. Glass ampules (26084-20), plastic ampules (260845-00), 100mL bottles (26084-42) of m-ColiBlue24 and agar plates (28052-15), afford maximum shelf life and ease of use for pre-measured and prepared media. Enumerate total coliforms and E. coli on one petri dish. Differentiate colonies easily - red and blue indicate total coliforms and blue specifies E. coli. Get optimal recovery of stressed and injured organisms.

No other pad is compatible with the product and the use of any pad other than those recommended below may result in unreliable, inaccurate data. m-ColiBlue24 should only be used with pads manufactured by the following companies/brands: Pall Corporation (e.g., Sterile Petri Dishes with Absorbant Pads), Private Label/OEM products manufactured by Pall, Gelman, Sartorius (cellulose and glass fiber pads), Sartorius BioSart products with the integrated Petri dish, and Millipore EO sterilized pads.

With ready-to-use m-ColiBlue24 Broth, you can:

Read and confirm results in 24 hours

Get superior sensitivity - 1 CFU/100mL

Conduct emergency testing and routine monitoring

Eliminate special equipment and ultraviolet lamp

Minimize background growth of non-coliform bacteria

Technische Daten

Material: Glass
Volumen / Packungsgröße: 100 mL