



Pipet, serological, sterile, polystyrene w/ plug, 10 mL, 500/pk

Artikel-Nr.: 2092868

EUR Preis: Kontakt

Anrufen für Liefertermin

Disposable polystyrene with plug in mouth end. Calibrated "to deliver/blow out".

Sterile, disposable plastic pipet with plug. Bulk packaged in 25/pack Color coded. Black ascending/descending graduations. Calibrated "To Deliver". Accuracy $\pm 2\%$.

Technische Daten

Graduations: 0.10 mL

Kapazität: 10 mL

Volumen / Packungsgröße: 500 /pk