



Be Right™



BEC/EDI BEC/EDI Basic Electrochem. Cell

Artikel-Nr.:

X51V001

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Basic Electrochem cell for use with EDI101 rotating disc electrode, includes XM110 Pt counter electrode

The BEC/EDI is a Basic Electrochemical Cell designed to support an EDI101 Rotating Disc Electrode. Its body is made of pyrex with a double jacketed device allowing thermostatic experiments to be performed. This cell can also be used with any type of solid electrode. The setup includes a calomel reference electrode and a platinum wire counter electrode.