



Test kit, molybdate, model MO-2

Artikel-Nr.: 1419301
EUR Preis: Kontakt
Verfügbar

Gefahrstoffe



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method

Effectively monitors corrosion inhibitors

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.2
	1
Koffertyp:	F
Messbereich:	0 - 8 mg/L Mo
Messbereich (2):	0 - 40 mg/L Mo
Messbereich (3):	NA
Methode:	Color disc/ Mercaptoacetic acid
Modell:	MO-2
Parameter:	Molybdate - As Molybdenum
Plattform:	Color Disc
Stromversorgung:	NA
Test-Kit Modell:	MO-2
Versandgewicht:	3