



Be Right™



DEHA Color Disc Test Kit, Model DH-1

Artikel-Nr.: 2168200

VERALTETER ARTIKEL
Dieser Artikel ist nicht mehr verfügbar.

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. The Model DH-1 test kit is a boiler/cooling water kit utilizing the color disc method. N, N-Diethylhydroxylamine (DEHA) replaces hydrazine and its derivatives as a oxygen scavenger in steam generation systems. Contains 50 tests.

Simple, fast, color disc method

Measures N, N-Diethylhydroxylamine (DEHA) directly in three different concentration ranges

To measure other oxygen scavenger compounds, use Prod. No. 2349500, Oxygen Scavenger Test Kit

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

| | |
|------------------------|-------------------------------|
| Anzahl Bestimmungen: | 100 |
| Kleinste Schrittweite: | 1 µg/L |
| | 5 µg/L |
| | 22.5 µg/L |
| Koffertyp: | R |
| Messbereich: | 1 - 65 µg/L DEHA |
| Messbereich (2): | 5 - 375 µg/L DEHA |
| Messbereich (3): | 22,5 - 1,700 µg/L DEHA |
| Methode: | Color disc/ Iron reduction |
| Parameter: | DEHA |
| Plattform: | Color Disc |
| Test-Kit Modell: | DH-1 |
| Versandgewicht: | 9 |