



Be Right™



Test-Kit Kohlendioxid CA-23 1,25-25/5-40/5-100

Artikel-Nr.: 143601

EUR Preis: Kontakt

Verfügbar

Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as CO₂, with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Technische Daten

Anzahl Bestimmungen:	200
Fußnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kleinste Schrittweite:	1.25 2 5
Koffertyp:	D
Messbereich:	1.25 to 25 mg/L CO ₂
Messbereich (2):	2 to 40 mg/L CO ₂
Messbereich (3):	5 to 100 mg/L CO ₂
Methode:	Drop count titration/ Sodium hydroxide Phenolphthalein
Parameter:	Carbon Dioxide - As CO ₂
Plattform:	Tropfenzähltest
Test-Kit Modell:	CA-23
Versandgewicht:	1