



Be Right™



Test kit cyanide, model CYN-3

Artikel-Nr.: 201002

EUR Preis: Kontakt

Anrufen für Liefertermin

Gefahrstoffe



The main sources of cyanide, an extremely poisonous substance, are metal-cleaning and electroplating baths, gas scrubber, gas works, coke ovens, and various chemical processes. It can be removed easily from waste effluents through alkaline or neutral chlorination procedures. This kit provides a simple, low-cost method using a stepped color comparator.

Simple, low-cost method using a stepped color comparator

Achieve fast results with an easy-to-carry test kit

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.01
Koffertyp:	R
Messbereich:	0 - 0.3 mg/L CN-
Messbereich (2):	NA
Messbereich (3):	NA
Methode:	Color disc/ Pyridine-Pyrazolone
Modell:	CYN-3
Parameter:	Cyanide - As free CN-
Plattform:	Color Disc
Stromversorgung:	NA
Test-Kit Modell:	CYN-3
Versandgewicht:	2.5