



Be Right™



Chlor gesamt Test Kit CN-21P, 10 - 200 mg/L, 100 Stück

Artikel-Nr.: 2444400
EUR Preis: Kontakt
Verfügbar

Methode / Chemie: Tropfenzähltitration / Thiosulfat

Chlorine Test Kits are manufactured by Hach Company to ensure compliance with stringent quality control standards. Using the Drop Count Titration method, the CN-21P Chlorine Test Kit is used to determine high concentrations of chlorine in sanitizing applications. 10-200 mg/L; contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Technische Daten

Anzahl Bestimmungen:	100
Fußnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kleinste Schrittweite:	10
Koffertyp:	B
Messbereich:	-
Methode:	Drop count titration/ Thiosulfate
Parameter:	Chlorine, mid and high range - As total Cl ₂
Plattform:	Tropfenzähltest
Test-Kit Modell:	CN-21P
Versandgewicht:	1