



Be Right™



## Dissolved oxygen test kit, model OX-2P

Artikel-Nr.: 146900

EUR Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Low dissolved oxygen levels in surface water systems usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This test kit utilizes drop count titration and contains 100 tests.

Method / Chemistry: Drop count t

**Utilizes drop count titration - a simple, economical titration method**

---

### Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.2
Koffertyp:	D
Messbereich:	- O <sub>2</sub>
Messbereich (2):	1 - 20 mg/L O <sub>2</sub>
Methode:	Drop count titration/ Modified Winkler
Parameter:	Dissolved Oxygen - As DO
Plattform:	Tropfenzähltest
Test-Kit Modell:	OX-2P
Versandgewicht:	3