



Be Right™



Beaker, polypropylene, low form, 400 mL, 6/pk

Artikel-Nr.: 108078

EUR Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Technische Daten

Durchmesser:	99.1 mm with spout
Einheiten:	6/pk
Gehäusematerial:	Polypropylene
Kapazität:	400 mL
Material:	Polypropylene, Low Form,