



**Be Right™**



## **Pipet tips, Tensette pipette 1970010, 1.0-10.0 mL, non-sterile, pk/250**

**Artikel-Nr.:** 2199725

**EUR Preis:** Kontakt

Lieferung innerhalb von 2 Wochen

Perform accurate, economical pipetting with Hach's TenSette pipets. Single-channel, air-displacement pipets come in two models: 0.1-1.0 mL volumes (0.1 mL increments) or 1.0-10.0 mL (1.0 mL increments). Volume selection with a simple twist of the turret.

---

### **Technische Daten**

Volumen / Packungsgröße: 250 /pk