



Be Right™



Bag, sterile, WhirlPak, 207 mL, 500/pk

Artikel-Nr.: 2233100
EUR Preis: Kontakt
Lieferung innerhalb von 2 Wochen

Approved by U.S. Environmental Protection Agency (USEPA). Disposable, transparent polyethylene bags. Gas-sterilized, pull-tab openings. Can be marked with waterproof felt pen or marker. Sure "whirl" closure. Nasco.

Technische Daten

Kapazität: 207 mL
Volumen / Packungsgröße: 500