




Be Right™



Test tube/ MPN rack, polyethylene, 20 mm (O.D.), 40 holes

Artikel-Nr.: 2497902

EUR Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Twenty to ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Supports MPN tubes. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Technische Daten

Durchmesser:	20 mm
Küvetten Anordnung:	4 x 10
Küvetten-Einsätze:	40
Volumen / Packungsgröße:	1 each