



Be Right™



Bottle, wash, diamond alert, 500 mL, 6/pk

Artikel-Nr.: 2613049

EUR Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

Technische Daten

Kapazität: 500 mL

Material: LDPE

Volumen / Packungsgröße: 6 St.